

Wundnadelabmessungen

Dimensions of surgical needles

Dimensions des aiguilles à sutures




Dimensiones de las agujas de sutura



Empfohlen für Haut und Muskeln


Recommended for Skin and Muscles


Recommandé pour la peau et les muscles

Recomendado para cutáneas y musculares

Typ type tipo	Bezeichnung Description descripción	Größe size taille tamaño	Abmessung in mm size in mm dimension en mm tamaño en mm
A 	201	000	1,60 x 102
	202	00	1,60 x 93
	702	0	1,60 x 86
B 		1	1,40 x 80
	203	2	1,40 x 73
	204	3	1,30 x 68
	303	4	1,30 x 63
	304	5	1,20 x 58
	403	6	1,10 x 54
	404	7	1,10 x 50
	704	8	1,00 x 46
	804	9	1,00 x 42
	904	10	0,85 x 38
G 		11	0,85 x 34
	211	12	0,75 x 30
	212	13	0,75 x 26
	311	14	0,75 x 23
	312	15	0,65 x 20
	411	16	0,65 x 18
	412		
	712		
812			
912			

Typ type tipo	Bezeichnung Description descripción	Größe size taille tamaño	Abmessung in mm size in mm dimension en mm tamaño en mm
BE 	203	1	1,30 x 80
	204	2	1,30 x 73
	704	3	1,10 x 68
		4	1,10 x 63
GE 	211	5	1,10 x 58
	212	6	1,00 x 54
	712	7	1,00 x 52
		8	1,00 x 50
		9	0,85 x 46
		10	0,85 x 44
		11	0,85 x 42
		12	0,75 x 38
		13	0,75 x 36
		14	0,75 x 34
		15	0,75 x 32
		16	0,65 x 30
		17	0,65 x 28
		18	0,65 x 26
		19	0,65 x 24,5
		20	0,65 x 23
		21	0,55 x 20
		22	0,55 x 17,5

Typ type tipo	Bezeichnung Description descripción	Größe size taille tamaño	Abmessung in mm size in mm dimension en mm tamaño en mm
C 	205	00	1,40 x 73
	206	0	1,30 x 68
	706	1	1,30 x 63
		2	1,20 x 58
		3	1,10 x 50
		4	1,00 x 46
		5	1,00 x 42
		6	0,85 x 38
		7	0,85 x 34
	8	0,75 x 30	


Typ type tipo	Bezeichnung Description descripción	Größe size taille tamaño	Abmessung in mm size in mm dimension en mm tamaño en mm
L 	222	00	1,40 x 73
	322	0	1,30 x 63
		1	1,30 x 55
		2	1,10 x 50
		3	1,00 x 42
		4	0,85 x 34
		5	0,75 x 30
		6	0,75 x 26
	7	0,65 x 23	
	8	0,65 x 20	


Empfohlen für Haut und Muskeln


Recommended for Skin and Muscles


Recommandé pour la peau et les muscles


Recomendado para cutáneas y musculares


Typ type tipo	Bezeichnung Description descripción	Größe size taille tamaño	Abmessung in mm size in mm dimension en mm tamaño en mm
Ga 	213	000	2,00 x 90
	214	00	2,00 x 80
	313	0	1,80 x 72
	314	1	1,60 x 66
	413	2	1,60 x 60
	414	3	1,40 x 54
	714	4	1,40 x 46
	814	5	1,30 x 38
	914	6	1,10 x 32
		7	1,00 x 27
		8	1,00 x 23
	9	0,85 x 18	

Typ type tipo	Bezeichnung Description descripción	Größe size taille tamaño	Abmessung in mm size in mm dimension en mm tamaño en mm
H 	218	00	1,40 x 73
	318	0	1,30 x 63
		1	1,10 x 54
		2	1,00 x 42
		3	0,85 x 34
		4	0,75 x 26
	5	0,55 x 22	

Typ type tipo	Bezeichnung Description descripción	Größe size taille tamaño	Abmessung in mm size in mm dimension en mm tamaño en mm
König 	220	1	1,30 x 55
		2	1,30 x 50
		3	1,30 x 45
		4	1,30 x 41
		5	1,10 x 37
		6	1,10 x 34
		7	1,10 x 32
		8	0,75 x 30
		9	0,75 x 26
		10	0,65 x 23
		11	0,65 x 20
		12	0,65 x 18

Typ type tipo	Bezeichnung Description descripción	Größe size taille tamaño	Abmessung in mm size in mm dimension en mm tamaño en mm
Keith 	209	1	0,75 x 75
	210	2	0,75 x 70
	309	3	0,75 x 65
	310	4	0,75 x 60
	710	5	0,75 x 55
	810	6	0,75 x 50

Typ type tipo	Bezeichnung Description descripción	Größe size taille tamaño	Abmessung in mm size in mm dimension en mm tamaño en mm
E 	319	000	1,10 x 72
	320	00	1,10 x 63
	419	0	1,00 x 54
	420	1	1,00 x 46
	820	2	0,85 x 38
	920	3	0,75 x 30
		4	0,65 x 24
		5	0,55 x 22
		6	0,55 x 17,5

Typ type tipo	Bezeichnung Description descripción	Größe size taille tamaño	Abmessung in mm size in mm dimension en mm tamaño en mm
Mayo 	323	1	1,30 x 51
	324	2	1,30 x 48
	423	3	1,30 x 45
	424	4	1,30 x 42